



CALCUTTA SCHOOL OF TROPICAL MEDICINE

GOVERNMENT OF WEST BENGAL

108 Chittaranjan Avenue, Kolkata - 700073

Phone : 033 2212 3695/96/97 Fax : 033 2212 3698

Website : stmkolkata.org

Memo No. STM/STORE/QUOTE/92/2020-21 **QUOTATION NOTICE** Date : 08.01.2021

Sealed QUOTATIONS are invited from the Reputed Company/Direct Importers/Authorized distributor firms for supply of items as mentioned below at the Department of **Virology Unit, School of Tropical Medicine, Kolkata.**

The QUOTATION is to be submitted in the drop box which will be kept at the office chamber of the undersigned on all working days from **09.01.2021.**

The last date of submission of QUOTATION is **15.01.2021 up to 1 pm** and the same will be opened on the same date at 2 pm in presence of interested bidder, if any. The bidders will submit the QUOTATION in their own official pad along with the following essential documents :

- Valid Trade License / Enlistment
- PAN card of the Bidder /Bidder company
- GST Registration certificate
- IT Return of any two of last three financial year
- Valid P.Tax Certificate

Physical verification of samples and or physical demonstration may be done at the discretion of the hospital authority by a team of expert to adjudge the applicability or suitability of the product or products in the functional requirement of the patients since the items are life savings and the technical evaluation will be done accordingly in close observation of the product or products, technical brochures along with Test certificate produced and physical sampling thereof.

Item Name

UV Sterilization Box
(Specification enclosed in Annexure I)

Subha
08.01.2021

Medical Superintendent cum Vice Principal
School of Tropical Medicine, Kolkata

Medical Superintendent cum Vice-Principal
Central Hospital for Tropical Diseases
School of Tropical Medicine
Kolkata

Annexure I

Specifications for UV Sterilization Box

Easy to clean and easy to install

UV light

LED lights for alarm

Automatic UV switch off

Auto UV sterilization cycle

Metallic interior and exterior furniture in stainless steel

Resistant to corrosive chemicals or cleaning reagents

UV Wavelength: 270 - 290 nm or higher

Radiation Power: UVC 360 microwatts/CM² or compatible

Radiation dose per minute: 4.5-5 KJ

Sterilization Method: UV-C Light

Timer: UV timer to be present

Auto UV sterilization cycle:

Sterilization standard time: 95-99% in 2-30 mins

Indicator: LED lights for indications and alarms

Automatic UV-C light switch off when door is open

Power supply: 220-250 Volts/50 Hz.